

Ruby Mobile Shower System

Ref: 3017.00UK

Technical Specifications



Dimensions

- Overall height: 106cm
- Overall width: 94cm
- Overall depth: 58cm

Description

- Steel frame, white colour, anti-corrosion coated
- Covers in ABS PMMA
- x6 Ø 100mm directional wheels, 2 with brakes
- 2 stainless steel handrails
- Fresh water tank in stainless steel, capacity 23L
- Waste water tank in stainless steel, capacity 23L
- 24V shower pump, 11 litres/min

- 24V drainage pump, 14 litres/min
- Simultaneous pump function
- Temperature controller with double display (real temperature & set temperature)
- Red warning light (water too hot)
- Resistance 230V 2700W
- Integrated differential circuit breaker 20A 30mA
- Electric components inside IPX5 waterproof box
- IP68 waterproof light buttons
- Handheld shower with flow adjustable trigger
- Concealed disinfectant tank
- Silicone autoclave piping

- One bed protection sheet, four clips & one towel included

Conformity

- 93/42/CEE Medical Device - Class 1
- IPX4 - Class 1 - CEI 60601
- CE marking

Accessories

- 3017.10 Shower sheet for bed - 90cm
- 3017.11 Shower sheet for bed - 120cm
- 3017.20 White bath towel
- Ruby tap filling adapter

To order or enquire about **Ruby**, please use the contact details below:

Reval Continuing Care Ltd
Unit 17 Brynmenyn Business Centre
St Theodore's Way
Brynmenyn
Bridgend CF32 9TZ

Tel: +44 (0)1656 721 001
Fax: +44 (0)1656 720 006
Email: sales@revalcc.com
Web: www.revalcc.co.uk

REVAL
SUSTAINABLE WELL-BEING

Innovative equipment designed with people in mind

WRAS APPROVED

DO8 APPROVED

BHTA
Medical & Healthcare Association

TMV
SCHEMES

ISO 9001:2008
REGISTERED FIRM



New Ruby Mobile Shower System

*Bringing the bathroom
to the bedroom*



The NEW *Ruby* mobile shower system has been designed to shower patients who cannot access the bathroom independently or with assistance.

Whether incapacitated due to a critical condition, long term illness, or obesity, *Ruby* provides unique and uncompromising hygiene levels within the safety and comfort of the patients' own bed.

A lightweight and waterproof sheet is positioned underneath the patient, which then fastens to the bed's side rails, head, and foot boards. The waterproof sheet accommodates the patient, whilst ensuring full water containment from the shower hose, and can be laundered at 70°C for re-use.

Alternatively, a low cost replacement sheet can be used as a patient specific item, or simply disposed of. There is an integrated water vacuum which transfers waste water from the sheet into a dedicated *waste water tank* on the appliance, which can be drained after use.

Step One



Firstly, fill *Ruby's* warm water container.

Step Two



Position *Ruby* at the end of the bed and select temperature. Press the heater button and wait for *Ruby* to reach desired temperature (this takes approximately ten minutes from cold).

Step Three



Lay one half of the purpose designed, double-width towel over the mattress, and allow the other half to rest over the side rails.

Step Four



Next, fit the waterproof sheet under the patient. Position the sheet as evenly as possible using the fasteners provided, and fix to the side rails and bed ends to ensure water containment.

The **Ruby** is compact and highly manoeuvrable, and is supplied with a purpose designed, double-width drying towel and waterproof shower sheet. All components are fully integrated into the appliance making it safe and user friendly. The **Ruby** incorporates anti-scald protection by thermostatically controlled water, which can be pre-set and heated to suit the patient's specific needs. The control panel has a simple interface, whilst the shower hose has a trigger, making the water pressure easy to control and comfortable for the patient or care worker. The **Ruby** holds up to 23 Litres, which is sufficient to shower a patient.

All surfaces of the **Ruby** are smooth, making it easy to clean and disinfect, using its built-in hand operated disinfection applicator, located discretely at the front of the appliance. Reval produces its own disinfection cleaning chemical for use when cleaning the **Ruby**. All hoses are easy to remove and can be thermally sterilised and replaced in minimal time.



Step Five



Unclip the shower hose in readiness to use. Then unclip the water vacuum hose, and position it at the lowest point of the mattress where water pools fastest. Remove the patient's clothing and check the waters output temperature prior to application.

Step Six



When the water temperature has reached the desired pre-set temperature it can then be applied to the patient.

Step Seven



With the shower procedure complete the patient can be dried and the sheet removed and stored. The patient can then be dressed.

Step Eight



With the hygiene procedure complete **Ruby** can then be emptied of waste water, and cleaned in readiness for the next use*.



* If **Ruby** is to be stored away then we advise the warm water container be drained down. This is done in the same way as the waste water, using its dedicated drainage hose situated at the back of the appliance.